

# DNAnexus Command Line Interface

*Stephen Hartley, PhD*

*2019-12-02*

# Contents

<b>1</b>	<b>Setup</b>	<b>2</b>
<b>2</b>	<b>The DNANexus Web Interface:</b>	<b>5</b>
<b>3</b>	<b>A Simple Example:</b>	<b>6</b>
<b>4</b>	<b>The App Builder Wizard:</b>	<b>7</b>
<b>5</b>	<b>Modifying the App:</b>	<b>8</b>
<b>6</b>	<b>Building the App:</b>	<b>12</b>
<b>7</b>	<b>Uploading Input Files</b>	<b>13</b>
<b>8</b>	<b>Running the software</b>	<b>14</b>
<b>9</b>	<b>Downloading results</b>	<b>15</b>
<b>10</b>	<b>Other useful commands:</b>	<b>17</b>
<b>11</b>	<b>Legal</b>	<b>18</b>

# Chapter 1

## Setup

First we have to load the DNAnexus command line interface utility. Both CCAD and Biowulf have it installed as a module.

```
#CCAD:
module load dx
#BIOWULF:
#module load DNAnexus
```

Next we have to login, select a project, and

```
#Log in (only have to do this the first time)
#dx login
```

```
#See what projects you have access to:
dx find projects
```

```
#Select the project you want from the above list:
dx select project-FbvPXyQ0pgP1PQbkJG1vfQpz
```

```
## project-FbvPXyQ0pgP1PQbkJG1vfQpz : Steve@NCI (ADMINISTER)
## project-FXfKVQj99yPxP6YqB7ykfBZ8 : AzureWest_Project (CONTRIBUTE)
## project-Ff54zy09462pQB99BpVZv7kG : steveTestProject (ADMINISTER)
## project-FYf1F7j9BppKqvGPK5GvqBzz : Genetic Susceptibility to Subsequent Neoplasms (CONTRIBUTE)
## Selected project project-FbvPXyQ0pgP1PQbkJG1vfQpz
```

Once dx is loaded and a project is selected, you can use many standard linux commands on the project by prefacing the command with “dx”

```
dx ls
dx cd test
dx cd ..

dx mkdir test2
dx rmdir test2
```

```
## test/
## findNi
```

For a full list of available commands:

```
dx help all
```

```
## dx all: All commands
##
## Commands:
```

```

##
## add                Add one or more items to a list
## add developers     Add developers for an app
## add member         Grant a user membership to an org
## add stage          Add a stage to a workflow
## add users          Add authorized users for an app
## add_types          Add types to a data object
## api                Call an API method
## build              Upload and build a new applet/app, or a workflow
## build_asset        Build an asset bundle
## cat                Print file(s) to stdout
## cd                 Change the current working directory
## clearenv           Clears all environment variables set by dx
## close              Close data object(s)
## cp                Copy objects and/or folders between different projects
## describe           Describe a remote object
## download           Download file(s)
## env                Print all environment variables in use
## exit               Exit out of the interactive shell
## find               Search functionality over various DNAnexus entities
## find analyses      List analyses in the current project
## find apps          List available apps
## find data          List data objects in the current project
## find executions    List executions (jobs and analyses) in the current project
## find globalworkflows List available global workflows
## find jobs          List jobs in the current project
## find org           List entities within a specific org.           "dx find org
##                   members" lists members in the
##                   specified org           "dx find org
##                   projects" lists projects billed to the
##                   specified org           "dx find org
##                   apps" lists apps billed to the
##                   specified org Please execute "dx find
##                   org -h" for more information.
## find org apps      List apps billed to the specified org
## find org members    List members in the specified org
## find org projects   List projects billed to the specified org
## find orgs           List orgs
## find projects       List projects
## generate_batch_inputs Generate a batch plan (one or more TSV files) for batch
##                   execution
## get                Download records, apps, applets, workflows, and files
## get_details         Get details of a data object
## head               Print part of a file
## help               Display help messages and dx commands by category
## install            Install an app
## invite             Invite another user to a project or make it public
## list               Print the members of a list
## list database       List entities associated with a specific database. For
##                   example,           "dx list database
##                   files" lists database files associated
##                   with a specific database.
##                   Please execute "dx list database -h"
##                   for more information.
## list database files List files associated with a specific database
## list developers     List developers for an app
## list stages         List the stages in a workflow

```

## list users	List authorized users for an app
## login	Log in (interactively or with an existing API token)
## logout	Log out and remove credentials
## ls	List folders and/or objects in a folder
## make_download_url	Create a file download link for sharing
## mkdir	Create a new folder
## mv	Move or rename objects and/or folders inside a project
## new	Create a new project or data object
## new org	Create new org
## new project	Create a new project
## new record	Create a new record
## new user	Create a new user account
## new workflow	Create a new workflow
## publish	Publish an app or a global workflow
## pwd	Print current working directory
## remove	Remove one or more items to a list
## remove developers	Remove developers for an app
## remove member	Revoke the org membership of a user
## remove stage	Remove a stage from a workflow
## remove users	Remove authorized users for an app
## remove_types	Remove types from a data object
## rename	Rename a project or data object
## rm	Remove data objects and folders
## rmdir	Remove a folder
## rmproject	Delete a project
## run	Run an applet, app, or workflow
## select	List and select a project to switch to
## set_details	Set details on a data object
## set_properties	Set properties of a project, data object, or execution
## set_visibility	Set visibility on a data object
## setenv	Sets environment variables for the session
## ssh	Connect to a running job via SSH
## ssh_config	Configure SSH keys for your DNAnexus account
## tag	Tag a project, data object, or execution
## terminate	Terminate jobs or analyses
## tree	List folders and objects in a tree
## uninstall	Uninstall an app
## uninvite	Revoke others' permissions on a project you administer
## unset_properties	Unset properties of a project, data object, or execution
## untag	Untag a project, data object, or execution
## update	Update certain types of metadata
## update member	Update the membership of a user in an org
## update org	Update information about an org
## update project	Updates a specified project with the specified options
## update stage	Update the metadata for a stage in a workflow
## update workflow	Update the metadata for a workflow
## upgrade	Upgrade dx-toolkit (the DNAnexus SDK and this program)
## upload	Upload file(s) or directory
## wait	Wait for data object(s) to close or job(s) to finish
## watch	Watch logs of a job and its subjobs
## whoami	Print the username of the current user

## Chapter 2

# The DNANexus Web Interface:




Figure 2.1: A figure caption.

## Chapter 3

# A Simple Example:

We are going to do a very simple example in which we upload some data, run a program on it, and then download the result.

Say we want to extract all the dictionary words that begin with “ni”.

## Chapter 4

# The App Builder Wizard:

The first step is to create a new App. You can generate the files from a template, write them yourself, or you can use the dx app creation wizard. For your first project, I recommend the last option.

Since the wizard is interactive, I have already performed this step and the example directory contains a copy of the output:

```
cp -R ../examples/findNi_original ./ex/findNi
```



## Chapter 5

# Modifying the App:

The app produced by the wizard is very minimal, it doesn't actually DO anything.

```
ls ./ex/findNi/
```

```
## dxapp.json
## findNi_original
## Readme.developer.md
## Readme.md
## resources
## src
## test
```

```
ls ./ex/findNi/src
```

```
## findNi.sh
```

The main files that you will need to modify are the dxapp.json file:

```
cat ./ex/findNi/dxapp.json
```

```
## {
##   "name": "findNi",
##   "title": "Find Ni",
##   "summary": "Finds words that begin with ni",
##   "dxapi": "1.0.0",
##   "version": "0.0.1",
##   "inputSpec": [
##     {
##       "name": "wordsFile",
##       "label": "File of all words",
##       "class": "file",
##       "optional": false,
##       "patterns": [
##         "*"
##       ],
##       "help": ""
##     }
##   ],
##   "outputSpec": [
##     {
##       "name": "outfile",
##       "label": "Output file of words",
##       "class": "file",
##       "patterns": [
```

```

##      "*"
##    ],
##    "help": ""
##  }
## ],
## "runSpec": {
##   "timeoutPolicy": {
##     "*": {
##       "hours": 48
##     }
##   },
##   "interpreter": "bash",
##   "release": "16.04",
##   "distribution": "Ubuntu",
##   "file": "src/findNi.sh",
##   "execDepends": [
##     {"name": "openjdk-8-jre-headless",
##      "package_manager": "apt"},
##     {"name": "tabix",
##      "package_manager": "apt"}
##   ]
## },
## "regionalOptions": {
##   "aws:us-east-1": {
##     "systemRequirements": {
##       "*": {
##         "instanceType": "mem1_ssd1_x4"
##       }
##     }
##   }
## }
## }
## }

```

and also the bash script in the src directory:

```
cat ./ex/findNi/src/findNi.sh
```

```

## #!/bin/bash
## # helloWorldTestApp 0.0.2
##
## main() {
##
##   #First, let's print out the input file, to make sure the inputs
##   # got set properly:
##   echo "Value of wordsFile: '$wordsFile'"
##   echo "Value of wordsFile_name: $wordsFile_name"
##
##   #Next let's download the words file from the cloud and into
##   # our virtual machine.
##   dx download "$wordsFile" -o $wordsFile_name
##
##   #This part is the program itself. It takes the input and grabs
##   # all the words that begin with Ni
##   cat $wordsFile_name | grep "^ni" > words.that.begin.with.ni.txt
##
##   #This part uploads the resultant file to the cloud
##   # If you don't do this, then all the files you made will be

```

```
##      # deleted when the virtual machine gets shut down.
##      #Note: dx upload --brief returns a file code that
##      # the next tool needs in order to register the output.
##      outfile_dx_code=$(dx upload words.that.begin.with.ni.txt --brief)
##
##      #Finally: this part registers the output and connects it to this job
##      # That way the output will be accessible from the job screen.
##      dx-jobutil-add-output outfile "$outfile_dx_code" --class=file
##
##
## }
```

I have already generated new versions of these files.

For the json file, I just added some dependencies just to show how one would do that. In this case I told the VM to install tabix and openJDK. Note: I don't actually need them in this example.

```
cp ../examples/files/dxapp.v2.json ./ex/findNi/dxapp.json
cat ./ex/findNi/dxapp.json
```

```
## {
##   "name": "findNi",
##   "title": "Find Ni",
##   "summary": "Finds words that begin with ni",
##   "dxapi": "1.0.0",
##   "version": "0.0.1",
##   "inputSpec": [
##     {
##       "name": "wordsFile",
##       "label": "File of all words",
##       "class": "file",
##       "optional": false,
##       "patterns": [
##         "*"
##       ],
##       "help": ""
##     }
##   ],
##   "outputSpec": [
##     {
##       "name": "outfile",
##       "label": "Output file of words",
##       "class": "file",
##       "patterns": [
##         "*"
##       ],
##       "help": ""
##     }
##   ],
##   "runSpec": {
##     "timeoutPolicy": {
##       "*": {
##         "hours": 48
##       }
##     },
##     "interpreter": "bash",
##     "release": "16.04",
##     "distribution": "Ubuntu",
##     "file": "src/findNi.sh",
```

```
##      "execDepends": [
##          {"name": "openjdk-8-jre-headless",
##            "package_manager": "apt"},
##
##          {"name": "tabix",
##            "package_manager": "apt"}
##      ]
##  },
##  "regionalOptions": {
##      "aws:us-east-1": {
##          "systemRequirements": {
##              "*": {
##                  "instanceType": "mem1_ssd1_x4"
##              }
##          }
##      }
##  }
##  }
## }
```

Next we write the script itself. The main function gets executed once the virtual machine is spun up and everything is installed and prepped. Note that it sets a bunch of environment variables that link to the input files.

The program below runs through the provided words file and returns a file containing all the words that begin with “ni”.

```
cp ../examples/files/findNi.v2.bash ./ex/findNi/src/findNi.sh
cat ./ex/findNi/src/findNi.sh
```

```
## #!/bin/bash
## # helloWorldTestApp 0.0.2
##
## main() {
##
##     #First, let's print out the input file, to make sure the inputs
##     # got set properly:
##     echo "Value of wordsFile: '$wordsFile'"
##     echo "Value of wordsFile_name: $wordsFile_name"
##
##     #Next let's download the words file from the cloud and into
##     # our virtual machine.
##     dx download "$wordsFile" -o $wordsFile_name
##
##     #This part is the program itself. It takes the input and grabs
##     # all the words that begin with Ni
##     cat $wordsFile_name | grep "^ni" > words.that.begin.with.ni.txt
##
##     #This part uploads the resultant file to the cloud
##     # If you don't do this, then all the files you made will be
##     # deleted when the virtual machine gets shut down.
##     #Note: dx upload --brief returns a file code that
##     # the next tool needs in order to register the output.
##     outfile_dx_code=$(dx upload words.that.begin.with.ni.txt --brief)
##
##     #Finally: this part registers the output and connects it to this job
##     # That way the output will be accessible from the job screen.
##     dx-jobutil-add-output outfile "$outfile_dx_code" --class=file
##
## }
```

## Chapter 6

# Building the App:

Now we have to “build” the app, which packages up the app and uploads it to your current project

```
cd ./ex/  
dx build -f findNi
```

```
## WARNING:dxpy:name "findNi" should be all lowercase  
## INFO:dxpy:Deleting applet(s) applet-FgXgj5Q0pgPF3X4v4ZQv0j80  
## {"id": "applet-FgXgjFj0pgP3yGXV15Q0bBY9"}
```

## Chapter 7

# Uploading Input Files

Before we can run our new app, first we need to upload an input file. For this we use the upload command:

```
dx upload ../examples/indata/linux.words.txt --brief --path /test/linux.words.txt
```

```
## file-FgXgjG00pgPBV2bJ8Z7ZzY7g
```

## Chapter 8

# Running the software

Before we can run our new app, first we need to upload an input file. For this we use the upload command:

```
dx run findNi -i wordsFile=/test/linux.words.txt --destination /test/
```

## Chapter 9

# Downloading results

We can now download the results using the dx download command:

```
dx download /test/words.that.begin.with.ni.txt
```

*#Let's see the first 50 words:*

```
cat words.that.begin.with.ni.txt | head -n50
```

```
## Error: path "/mnt/nfs/gigantor/ifs/Shared/hartleys/docs/DnaNexusTutorial/build
## /words.that.begin.with.ni.txt" already exists but -f/--overwrite was not set
## ni
## niacin
## niacinamide
## niacins
## niagara
## niagra
## niais
## niaiserie
## nialamide
## niata
## nib
## nibbana
## nibbed
## nibber
## nibbing
## nibble
## nibbled
## nibbler
## nibblers
## nibbles
## nibbling
## nibblingly
## nibby
## nibby-jibby
## nibelung
## niblic
## niblick
## niblicks
## niblike
## nibong
## nibs
## nibsome
## nibung
```



## nicad  
## nicads  
## nicaean  
## nicaragua  
## nicaraguan  
## nicaraguans  
## niccolic  
## niccoliferous  
## niccolite  
## niccolo  
## niccolous  
## nice  
## niceish  
## niceling  
## nicely  
## nicene  
## nice-Nellie

# Chapter 10

## Other useful commands:

You can use cat to quickly preview files:

```
dx cat /test/words.that.begin.with.ni.txt | head -n10
```

```
## ni
## niacin
## niacinamide
## niacins
## niagara
## niagra
## niais
## niaiserie
## nialamide
## niata
```

```
set -x
```

```
dx describe /test/words.that.begin.with.ni.txt
```

```
## + dx describe /test/words.that.begin.with.ni.txt
## Result 1:
## ID                file-FgXgbfj0gzgbqx5KJ97V6J7j
## Class             file
## Project           project-FbvPXyQ0pgP1PQbkJG1vfQpz
## Folder            /test
## Name              words.that.begin.with.ni.txt
## State             closed
## Visibility        visible
## Types             -
## Properties        -
## Tags              -
## Outgoing links    -
## Created           Mon Dec  2 15:58:03 2019
## Created by        hartleys
## via the job       job-FgXgbK80pgPJJ54XF2qK4zBp
## Last modified     Mon Dec  2 15:58:05 2019
## archivalState     "live"
## Media type        text/plain
## Size              10.33 KB
dx watch <jobid>
```

# Chapter 11

## Legal

This document and all relevant documentation is “United States Government Work” under the terms of the United States Copyright Act. It was written as part of the authors’ official duties for the United States Government and thus vArmyKnife cannot be copyrighted. This software is freely available to the public for use without a copyright notice. Restrictions cannot be placed on its present or future use.

Although all reasonable efforts have been taken to ensure the accuracy and reliability of the software, documentation, and data, the National Human Genome Research Institute (NHGRI), the National Cancer Institute (NCI) and the U.S. Government does not and cannot warrant the performance or results that may be obtained by using this software or data. NHGRI, NCI and the U.S. Government disclaims all warranties as to performance, merchantability or fitness for any particular purpose.

In any work or product derived from this material, proper attribution of the authors as the source of the software or data may be made, using “NCI Division of Cancer Epidemiology and Genetics, Human Genetics Program” as the citation.